

Club Warehouse Sizing Chart

When applying to limbs, measure to fit the narrowest part. When applying to other areas of the body measure the affected area

****Guide is for Medium Compression; for Low Compression select next size bandage up**

Measuring Guide for Application

Size	Application	Roll	Circumference Measurement**
AA	Infant Limbs	10m	< 6cm*
A	Infant Limbs	10m	7 – 11cm*
B	Small Limbs	10m	12 – 14cm*
C	Adult Limbs	10m	15 – 25cm*
D	Large Limbs	10m	26 – 35cm*
E	Adult Legs	10m	36 – 44cm*
F	Large Legs	10m	45 – 50cm*
G	Large Thighs	10m	51 – 60cm*
H	Extra Large Thighs	10m	61 – 70cm*
J	Small Trunk	10m	71 – 75cm*
K	Medium Trunk	10m	76 – 80cm*
K2	Medium Large Trunk	10m	81 – 84cm*
L	Large Trunk	10m	85 – 90cm*
M	Extra Large Trunk	10m	> 90cm*

